MR 3003 control

**Betriebsanleitung
Instruction Manual
Mode d'Emploi
Instrucciones de
Empleo**

DEUTSCH

D

ENGLISH

E

FRANCAISE

F

ESPAGNOL

ES

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- EN 61 010-1
 EN 61 010-1/A2: 1995

Thank you for buying a Heidolph Instruments product. This item has been designed, made and inspected in compliance with DIN EN ISO 61010 for long-term and flawless operation.

SUMMARY

SUMMARY	20
STANDARD HARDWARE & OPTIONS	21
GENERAL INFORMATION	22
SAFETY INFORMATION	22
SET-UP	23
1. machine set-up	23
2. electric hook-up	23
3. connect external temperature sensor Pt 1000 (in chromium-steel or glass sheath).....	23
4. connect external 2-station temperature sensor Pt 1000 (in chromium steel- or glass sheath).....	24
5. set-up using tank adapter	24
OPERATION AND CONTROLS	24
1. turn item ON	24
2. set stirring speed	25
3. heater control using temperature sensor (Set1 mode)	25
4. heater control using temperature sensor for 2 different stations (Set2 mode).....	27
5. heater control in absence of temperature sensor (Set3 mode)	28
6. residual heat indicator.....	28
7. interface installation	29
8. error messages and corrective action.....	30
INFORMATION ABOUT USING TEMPERATURE SENSORS	31
CLEANING & SERVICING	33
UNINSTALL, FORWARD & STORE	33
DISPOSAL.....	33
MALFUNCTION.....	33
SPECIFICATIONS.....	34
WARRANTY, LIABILITY AND COPYRIGHT.....	35
QUESTIONS / REPAIR WORK.....	35
CE-DECLARATION OF CONFORMITY.....	35

STANDARD HARDWARE & OPTIONS

Standard Hardware

PRODUCT	QUANTITY	P/N	
		230 V / 50/60 Hz	115 V / 60 Hz
MR 3003 control S or	1	504-33200-00	504-33200-01
MR 3003 control C or	1	504-33210-00	504-33210-01
MR 3003 control G	1	504-33220-00	504-33220-01
Instruction Manual	1	01-005-002-90	01-005-002-90
Power cord	1	14-007-003-81	14-007-003-89




Options

PRODUCT	P/N
Temperature sensor Pt 1000, chromium steel sheath	509-67910-00
Temperature sensor Pt 1000, glass sheath	509-67920-00
Temperature sensor Pt 1000, includes 2 measuring stations, chromium steel sheath	509-67930-00
Temperature sensor Pt 1000, includes 2 measuring stations, glass sheath	509-67940-00
1 liter tank - for oil up to 250 °C max.	504-93000-00
1 liter bath, PTFE-coated - for water up to 100 °C max.	504-93100-00
2 liter bath - for oil up to 250 °C max.	504-92000-00
2 liter bath, PTFE-coated - for water up to 100 °C max.	504-92100-00
4 liter bath - for oil up to 250 °C max.	504-91000-00
4 liter bath, PTFE-coated - for water up to 100 °C max.	504-91100-00
Adapter ring for 1 liter round-bottom flasks, non-coated	504-94000-00
Stand, 1.4305 (V2A) stainless steel	509-81000-00
Stand, 1.4305 (V2A) stainless steel (for 4 litre bath)	509-98000-00
Temperature sensor holder	509-67000-00
Magnetic paddle kit (3 ea.), cylindrical, PTFE-sheath	509-56000-00



GENERAL INFORMATION

-  Unpack your item carefully. Inspect for damage and report such damage or missing parts to your supplier right away.
-  Read your instruction manual carefully. Take time to save time while trying to work with your item. Make sure that every user has read and understood the instruction manual.
-  Please locate the instruction manual in a place easily accessible to every user.

IF ALL ELSE FAILS, READ THESE INSTRUCTIONS !

-  A so called EURO- plug (DIN 49441 CEE 7/II 10/ 16 A 250 V) is standard on all of the items. For the Continental US they feature a US-standard plug (NEMA pub.no.WDI.1961 ASA C 73.1.1961 page 8 15A 125V).
-  For using the item in a country with deviating outlet / plug systems, we recommend to use approved adapters or to have an electrician replace the standard plug with one suiting your needs.
-  As shipped, the item features a protective ground wire. When replacing the original plug, make sure to re-connect this protective ground wire in the new plug !

SAFETY INFORMATION

-  Please comply with all safety and accident-prevention regulations as in force for laboratory work !
-  Use extra care when working with flammable substances; refer to safety data sheets.
-  When connecting your item with your local power supply, please make sure your item is designed for your local supply voltage; go by data plate on the item.
-  Turn your power switch OFF whenever the item is not used, or before disconnecting the plug.
-  Repair work is limited to technicians approved by Heidolph Instruments.
-  Use extra care when working in the vicinity of flammable and explosive substances. Motors are of non-sparking type, the item itself however is not explosion-protected.
-  Please connect your magnetic stirrer with a protective-ground outlet only.
-  Caution! Hazard of serious burns when running the hot plate over 50 °C. Avoid contact with hot plate, bath and fluid.
-  Caution! When heating flammable substances, adjust hot plate nominal temperature to a lower level than the substance's flashpoint.
-  Your item needs a solid stand.

SET-UP


1. Set-up

- Please locate the item on a sturdy, horizontal surface.
- For reasons of safety, avoid arranging more hardware around the stirrer module.

2. Electric hook-up

front view, MR 3003 control

rear view, MR 3003 control

Do the following checks before connecting:

- Power switch (1) OFF (position "0"). (When ON, the toggle switch lights GREEN)
- Control knobs (2 + 3) in extreme left position (CCW)
- Connect item with mains by plugging the power cord in the plug connector (15) located in the item's rear panel and your outlet.

3. Connect external temperature sensor Pt 1000 (chromium steel- or glass sheath)

- Turn item OFF attempting to connect the temperature sensor.
- Connect sensor's bayonet-type plug with before flanged connector (16) in the item's rear panel

Please keep in mind that sensor needs an immersion depth of not less than 20 mm (mark) to avoid false reading.


For Pt 1000 wiring, please refer to diagram below.

Item does accept a Pt 1000 temperature sensor only (P/N 509-67910-00 or 509-67920-00); a Pt 100 won't work.

set-up includes stand and temperature sensor holder

3.1 Pt 1000 wiring

4. **Connect external 2-station temperature sensor Pt 1000 (chromium steel- or glass sheath)**
 - For set-up and controls refer to page 27.
5. **Set-up using bath**

If you are going to use a bath (refer to „options“ section) make sure, bottom ring of your adapter properly slips over the hot plate.

When using the item with a 4 liter bath, you may attach the temperature sensor to a special stay bar (18) (P/N 509-98000-00), designed for the 4 liter bath.


18

Set-up features 4 liter bath,
special stay bar and temperature sensor

OPERATION AND CONTROLS


1. Turn item ON

- Locate vessel or bath (to include stirring paddle) on the plate of your magnetic stirrer. Magnetic paddles are available from Heidolph Instruments as an extra option.

- *Prior to* turning the power switch (1) ON, rotate speed controls (2 + 3) to their extreme left position (CCW).
- Turning the power switch (1) ON (range I), a self-test is done; all LEDs (8, 11) light for a moment, and all bars in the digital display (7) appear, display shows **Set1**, **Set2** or **Set3**, followed by 3 beeps.
- When power switch is ON, stirrer motor is ready to run.
- **Act Temp** (11) LED lighting in **Set1** and **Set2** mode, display (7) will show actual sensor temperature, or, in **Set3** mode, hot plate temperature.

2. Set stirring speed


- Set required speed with the rotary speed control (2).
- To display actual speed, hit **Select** key (5) several times, until **Act Speed** LED (14) lights. Speed appears in display (7).

3. Heat medium, using temperature sensor (Set1 mode)

- Connect temperature sensor as described in para. 3 „Set-up“.

3.1 Set nominal temperature

- Hit **Select** key (5) several times, until **Set Temp** LED (12) lights.
- Thereafter hit **Set** key (6). **Set Temp** LED (12) starts flashing, display (7) shows temperature selected with rotary temperature control (3); at this point, you may alter temperature setting with rotary control (3).
 - ☞ Set temperature can only be changed after having interrupted the heating by pressing key **on/off** (4). LED **Heater** (8) does not light.
- To accept temperature, hold **Set** key (6) depressed until alteration is confirmed by a beep and „**Act Temp**“ appears in the display.
- Setting action may be discontinued by hitting **Set** key (6) once. This way, you return to the previous setting.
 - ☞ Attempting to lower nominal temperature by an amount exceeding ΔT set, will be reported as error code F6.
Remedy: rise ΔT -value until system temperature has dropped to an acceptable level.

3.2 Set safe temperature using ΔT (temperature difference)

Safe temperature is the total of temperature set plus ΔT (temperature cushion). Exceeding this safe temperature will disconnect your magnetic stirrer for reasons of safety.

- For adjusting ΔT hit **Select** key (5) several times, until **Set ΔT** (13) LED lights.
- Hit **Set** key(6), **Set ΔT** LED (13) starts flashing and display (7) shows temperature set with rotary temperature control (3).
 - ☞ Temperature difference ΔT can only be changed after having interrupted the heating by pressing key **on/off** (4). LED **Heater** (8) does not light.
- To accept temperature, hold **Set** key (6) depressed until alteration is confirmed by a beep and „Act Temp“ appears in the display.
- Setting action may be discontinued by hitting **Set** key (6) once. This way, you return to the previous setting.

3.3 Heating circuit ON / OFF

After having set speed, nominal temperature and ΔT , you may turn the heating circuit ON. At this point you may double-check set values by hitting the **Select** key (5) several times. This way, **Act Temp**; **Set Temp**; **Set ΔT** and **Act Speed** will appear in the display (7), one after the other.

- We recommend to select **Act Temp** display all the time.
- Hit the **ON/OFF** key (4) to turn the heating circuit ON.
- **Heater** LED (8) is going to light, i.e. system is heating.
- A yellow LED (9) confirms heating current ON.
- A soon as hot plate temperature exceeds 50°C with power switch or heating circuit at OFF, a RED LED (10) is starting to flash.
(also refer to para. 6, residual heat warning)
- Hit the **ON/OFF** key (4) to turn the heating circuit OFF.

☞ All of the parameters so set are saved by the built-in memory circuit and are available at re-start.

4. Heat medium, using 2-station temperature sensor (Set2 mode)

Connect 2-station temperature sensor as described in para. 3 „Set-up“.

For special applications, we recommend a 2-station temperature sensor (refer to picture). This sensor allows to limit max. temperature of the medium to a certain ΔT from nominal and avoids temperature overshoot in the medium close to the bath wall. For working with this temperature sensor please refer to the following details:


 Safety cut-off feature (refer to para. 3.2 „set safe temperature“) is not active.

 Heating time may rise at low ΔT .


 Medium temperature overshoot might occur at high ΔT due to the inherent thermal capacity of the heat carrier. If this happens, decrease ΔT .

- Immerse a temperature sensor with white handle (a) (= marked **medium**) in the medium bath.
- Immerse the second temperature sensor with black handle (B) (= marked **heat transfer medium or carrier**) in the heat transfer medium.

Use one more holder for this temperature sensor (P/N 509-67000-00).

set-up
including 2-station temperature sensor and heating bath

wiring
2-station temperature sensor

4.1 Select two-station temperature sensor control mode

Switch control mode from **Set1** to **Set2**.

- To this end turn power switch (1) OFF
- Depress and hold *both*, **Select** key (5) and **Set** key (6), turn power switch (1) ON.
- Release both keys
- Hit **Select** key (5) once, **Set2** appears in display (7).
- Hit **Set** key (6) to confirm **Set2** mode
- **Set2** mode ready is confirmed by 3 beeps and medium (substance) temperature appearing in the display (7). **Act Temp** LED (11) lights.

4.2 Setting 2-station temperature sensors

Set nominal temperature, ΔT and turn the heating circuit ON as described under para. 3.1, 3.2 and 3.3.

 Whenever you replace the 2-station temperature sensor with a „normal“ Pt 1000 sensor (chromium steel or glass sheath), change control mode to Set1 (refer to para. 4.1 „select 2-station temperature sensor control mode“)

5. Heating in absence of temperature sensor (Set3 mode)

No temperature sensor connected will make the item control hot plate temperature only. MR 3003 control is designed as to auto-detect this operational mode (**Set3** mode).

- When turning the power switch (1) ON, **Set3** mode detected appears in the display (7).
- After 3 beeps, display (7) shows the hot plate temperature with a prefix **H**; **Act Temp** LED(11) lights.
- For setting hot plate nominal temperature please refer to para. 3.1 „set nominal temperature“. Also refer to para. 3.3 „turn heating circuit ON/OFF“.

6. Residual heat indicator

To avoid burns when touching the hot plate, the MR 3003 control features a residual heat indicator. A flashing red LED (10) warns you that hot plate temperature still is above 50°C; this safety feature is available only with the power switch or heater switch in OFF position.

7. Interface connection

- Turn the item OFF before connecting an interface cable.
- We recommend to use an interface cable with SUB-D 9 connectors (male), use the female connector (17) located in the item's rear panel.


rear view, MR 3003 control

Do not use but the wires shown to avoid malfunction of your interface.

7.1 RS 232 interface

- This interface allows to define nominal temperature of your samples and nominal speed.
- This interface allows to send OUT actual temperature of your samples and actual speed.
- For wiring details of the SUB-D connector refer to drawing.
- A dot in the display is going to flash as long as the interface is active.

Caution: Do not use but the wires shown to avoid malfunction of your interface.

7.1.1 RS 232 interface commands:


trigger		query	
command	action	command	action
REM!	activate interface	TP1?	act temp Pt 1000 probe1
LOC!	deactivate interface	TP2?	act temp Pt 1000 probe2
H1!	heating circuit ON	TH1?	act temp Pt 100 hot plate1
H0!	heating circuit OFF	D?	act speed
M1!	select Set1 mode		
M2!	select Set2 mode		
D=xxxx!	Set speed		
TP=xxx!	Set temp probe (Set1 & Set2 only)		
DP=xxx!	Set ΔT (Set1 & Set2 only)		
TH=xxx!	Set hot plate temp (Set3 only)		

7.1.2 Interface parameters

- 9600 baud, no parity, 8 bit, 1 stop bit, no protocol

7.2 Analog interface 0 – 10V

- The analog interface allows to send OUT actual sample temperature.
- For wiring details of the SUB-D connector refer to drawing.
- Data sent OUT: 0 °C = 1 V; every 100 °C = 1 V

Caution:

Do not use but the wires shown to avoid malfunction of your interface.

8. Error messages and remedy

If one of the errors stipulated below is found, the microprocessor turns the MR 3003 control to OFF. TRIAC and relay are disconnected by the relay, the heating circuit enters stand-by mode, the error message appears in the display, and is confirmed by a beep.

Possible error messages:

- F2 Hot plate temperature more than 360 °C
 F3 Pt 100 in the hot plate are found beyond logical readings
 F4 Difference between Pt100 and hot plate exceeds 30 K
 F5 **Set1** mode: difference between both of the Pt 1000 exceeds 10 K
 F6 **Set1** mode: sample temperature overshoot by more than ΔT
 F7 Triac broken
 F8 **Set1**; **Set2** or **Set3** mode once detected has been changed during operation
Set1 mode, **Set2** mode: Pt 1000 disconnected or poor contact
Set3 mode: Pt 1000 connected
 F10 Processor failure

Turning the power switch ON/OFF will delete all error messages

If error messages F2 to F7 and F10 appear repeatedly, we recommend to have a skilled technician inspect your item.


INFORMATION ABOUT USING TEMPERATURE SENSORS

For more information about the function principle of temperature sensors please refer to items 3 & 4, para. „operation and controls“

1. How to work with a temperature sensor (Set1 mode)

The temperature covers several mission types.


1.1 Temperature sensor located in reactive substance, no heat transfer substance used

Benefits:

- Rapidly reaches an exact nominal temperature
- Adjustable safety margin

1.2 Temperature sensor in reactive substance, uses heat transfer substance


Benefits:

- Exact nominal temperature reached at ease
- Adjustable safety margin

Limitation:

- Allows a heat transfer substance temperature significantly above nominal temperature.

1.3 Temperature sensor in heat transfer substance

Benefits:


- Rapidly reaches nominal temperature
- Adjustable safety margin

Limitation:

- Substance temperature will stabilize on a random level below nominal, as a function of temperature in the heat transfer substance.

2. Use a 2-station temperature sensor (Set2 mode)

2.1 Temperature sensor “substance medium” in the reactive substance, temperature sensor “heat transfer medium” in the heat transfer substance

Benefits:

- Exact nominal temperature reached at ease
- No thermal load on thermolabile substances

Limitation:

- Reaches nominal value slowly

CLEANING & SERVICING

Cleaning: wipe housing clean with a damp cloth (add some sort of mild liquid soap).

Note

To avoid damage to the surface finish, avoid using chlorine bleach, chlorine-based detergents, abrasive substances, ammonia, rags or cleaning agents containing metal particles.

The item is maintenance-free. Repair work is limited to technicians approved or appointed by Heidolph Instruments. Please call your local Heidolph Instruments dealer or a Heidolph Instruments field representative (also refer to page 35).

UNINSTALL, FORWARD & STORE

Uninstall

Turn the item OFF and disconnect mains plug.

Forward & Store

1. We recommend to store the item and its components in its original box, or a similar container that offers adequate protection against damage in transit. Tape the box securely.
2. Store the item in a dry place.

Caution

Do not jolt or shake the item during transport.

DISPOSAL

For disposal, please comply with your local or national regulations.

Split by metal, plastic, etc.

Packing material to be treated as described above (material split).

MALFUNCTION

power switch dim:

- check power supply and power cord

won't stir:

- no magnetic paddle in stirring bath
- speed control knob in extreme left position (CCW)
- interface not connected as specified (refer to para. 7 „interface connection“)

won't heat:

- check error messages displayed (refer to para. 8 „error messages and remedy“)
- interface not connected as specified (refer to para. 7 „interface connection“)

F2 - F10 error displayed:

- refer to para. 8 „error messages and remedy“ under „operation & controls“

In case of a malfunction you cannot repair as described above, get in touch with your authorized Heidolph Instruments dealer.

SPECIFICATIONS

stirrer drive	capacitor motor, electronic speed control, thermal circuit breaker
speed range	0 – 1300 rpm
speed control	digital, 4-digit LED display, green
temperature control	electronic temperature control for hot plate and substance temperature, features microprocessor
heating power	600 watts
temperature range, hot plate	0 – 300 °C
hot plate made from	standard: Silumin (aluminum alloy), marked S chromium steel, marked C vitrified ceramics, marked G
temperature sensor in hot plate	2 x Pt 100
temperature range, substance	room temperature ⇒ approx. 250 °C
temperature range, substance	temperature sensor Pt 1000 (optional): 2 x Pt 1000, single-station temperature sensor Pt 1000 (optional): 2 x Pt 1000, 2 separate stations
control hysteresis, substance	< +/- 0,2 K
safety circuits, hot plate temperature	2 separate safety circuits featuring 2 contactors (triac and relays)
RS 232 serial interface	IN: nominal sample temperature and nominal speed; OUT: actual sample temperature and actual speed
analog interface (0 – 10 V)	OUT: actual sample temperature
supply voltage	230 volts / 50/60 Hz; others available on request
power dissipation	630 watts
weight	3.0 kg
dimensions (w x d x h)	240 mm x 155 mm x 120 mm
protection class	IP 30
ambient specifications	0 – 40 °C at 80 % rel. humidity
chemical stability	housing: encapsulating, polyamide coated; front panel: polyester

WARRANTY, LIABILITY & COPYRIGHT

Warranty

Heidolph Instruments warrants that the present product shall be free from defects in material (except wear parts) and workmanship for 3 years from the date shipped off the manufacturer's warehouse.

Transit damage is excluded from this warranty.

To obtain such warranty service, get in touch with Heidolph Instruments (phone: +49 – 9122 - 9920-68) or your local Heidolph Instruments Dealer. If defects in material or workmanship are found, your item will be repaired or replaced at no cost. Misuse, abuse, neglect or improper installation are not covered by this warranty promise. Alterations to the present warranty promise need Heidolph Instruments' consent in writing.

Exclusion Clause

Heidolph Instruments cannot be held liable for damage from improper use or misuse. Remedy for consequential damage is excluded.

Copyright

Copyright in pictures and wording of the present instruction manual is held by Heidolph Instruments.

QUESTIONS / REPAIR WORK

If any **aspect** of installation, operation or maintenance remains unanswered in the present manual, please get in touch with the following address:

For repair jobs please call Heidolph Instruments (phone: +49 – 9122 - 9920-68) or your local, authorized Heidolph Instruments dealer.

Note

You will receive clearance for sending your defective item to the following address:

Heidolph Instruments GmbH & Co. KG
Lab Equipment Sales
Walpersdorfer Str. 12
D-91126 Schwabach / Germany
Tel.: +49 – 9122 - 9920-68
Fax: +49 – 9122 - 9920-65
E-Mail: sales@heidolph.de

Safety Information

When shipping items for repair that may have been contaminated by hazardous substances, please:

- advise exact *substance*
- take adequate *protective* action not to harm our parts receiving and service personnel
- mark the pack IAW Hazardous Materials Act
-

CE-DECLARATION OF CONFORMITY

We herewith declare that the present product complies with the following standards and harmonized documents:

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© Heidolph Instruments GmbH & Co. KG